



The First Choice in
Vision Testing Systems

**Serving Eye Care
Professionals
Since 1992**

January 2011

Computer\LCD Requirements for Smart System Series Visual Acuity Systems

Computer: - PC or Apple, tower or mini tower or slim profile. ***If the user will connect more than 1 LCD screen, then a dual VGA graphics card is required and not supplied by M&S.*** Standard USB ports are required.

Laptop Computers: - can be used effectively for stand alone systems using the laptop display as the primary acuity display.

***** When attempting to use a laptop as a “central” cpu for operating dual monitors, it cannot be supported and it is NOT recommended.***

Processor: - Intel Celeron 1.8 Ghz or AMD Athalon 64 2.0Ghz or better

Operating System: - Windows XP with Service Pack 2 installed (SP2), Vista or Windows 7, Home or Professional

RAM Memory: - Windows® XP 512 meg, Windows® Vista and Windows® 7 is 2 Gigabyte

Hard Drive: - Minimum 40 Gigabyte

DVD Drive: - Required for video control via wireless remote control.

LCD: - Minimum 17”, 4 x 3 configuration or 16 x 9 configuration with 1280 x 1024 minimum resolution

Anything less than minimum required will not be supported by M&S.

Contact M&S Technologies if you have any questions.

www.mstech-eyes.com

Phone: (877) 225-6101, (847) 763-0500 Fax: (877) 473-9806, (847) 763-9170

M&S Technologies Inc., P.O. Box 1171 Park Ridge, IL 60068, Manufacturing Plant: 5557 W. Howard Street, Skokie, IL 60077