



## Smart System® | VR Headset

SMARTTracker™ 2 Advanced Eye Tracking Technology

### The Only VR Headset with Time-Tested Software

**M&S takes vision testing to new heights.** The programmable VR Headset with Surface Pro 7 Tablet is powered by our trusted M&S Smart System® software for fast and accurate vision testing. Our years of experience and service continue to help your clinic shine.



#### Reliable Results

Trust data accuracy with SMARTTracker™ 2, the most sensitive fixation monitoring technology available on the market today.



#### Speed and Efficiency

Our proven testing protocols reduce patient fatigue by providing QuickScreen results in 45 seconds and full threshold test results in under 3 minutes per eye.



#### Patient and Practice Friendly

Comfortable design with easy set up, audio and visual cues to your technician and patient, the ability to test without a wifi connection and results that can be immediately exported to your EMR.

### THE RIGHT TECHNOLOGY



847-763-0500  
mstech-eyes.com

The First Choice in Vision Testing Systems

A HILCO VISION COMPANY

# Innovative Technology, Proven Expertise

The Smart System® VR Headset is the newest addition to over 30-years of trusted quality and service from M&S Technologies. The programmable headset is lightweight, compact and portable, allowing your clinic to test in a variety of settings.

Combine the portability of the VR headset with the accuracy of our Smart System® software and you get the most versatile vision testing device on the market. Our team not only has decades of experience developing and refining vision testing software, but also provides the dependable, trained support needed to back it up.



## VR Headset Tests Included


### Visual Field

- SMARTracker™ 2 monitors and maintains patient fixation on central target
- 10-2, 24-2, 30-2
- 45 second QuickScreen;
- Full test results in under 3 minutes per eye
- Supra Threshold
- Ptosis
- Esterman

### Color Vision

### Stereo Testing


### Contrast Sensitivity



Visual Field Test Result

## Added Benefits

- Exportable results immediately available on Surface Pro 7 Tablet
- Progression analysis
- Visual field test pause / restart feature ensures proper eye fixation for greater accuracy
- Suitable for easy remote testing
- Lightweight, cushioned and compact design engineered for patient comfort
- Test in a fully illuminated room
- NO internet connection required
- English and Spanish audio instructions
- Live, immediate technical support team



SMARTracker™ 2  
*The most advanced eye tracking software on the market.*



ADA Compliant

**847-763-0500**  
[mstech-eyes.com](http://mstech-eyes.com)

The First Choice in Vision Testing Systems

©2022 M&S and Smart System are registered trademarks of the Hilsinger Company Parent LLC. All rights reserved. The Hilsinger Company Parent LLC holds US Patents 7,354,155; 7,926,948; 8,425,040; 8,167,429; 8,419,184; 8,550,631; 8,992,022; 9,433,347; 9,820,644; 10,244,938 and 10,182,713. Other Patents Pending.



A HILCOVISION COMPANY